

Product Datasheet

PAR-2 polyclonal antibody (orb1089101)

Description	PAR-2 polyclonal antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human PAR-2:37-100/397
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~40kDa
Uniprot ID	P55085
Dilution Range	IHC-P=1:100-500 □□ not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Expiration Date	12 months from date of receipt.