

Product Datasheet

Glypican 1 polyclonal antibody (orb1089123)

Description	Glypican 1 polyclonal antibody
Species/Host	Rabbit
Reactivity	Bovine, Canine, Equine, Mouse, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	KLH conjugated synthetic peptide derived from mouse Glypican 1 :251-350/557
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~62kDa
Uniprot ID	Q9QZF2
Dilution Range	Flow-Cyt=2ug/Test not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Expiration Date	12 months from date of receipt.