

Product Datasheet

Phospho-Estrogen Receptor alpha (Ser118) monoclonal antibody (orb1089146)

Description	Phospho-Estrogen Receptor alpha (Ser118) monoclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human ER alpha around the phosphorylation site of Ser118:QL(p-S)PF
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
MW	~66 kDa
Uniprot ID	P03372
Dilution Range	IHC-P=1:50-200 not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Expiration Date	12 months from date of receipt.