

Product Datasheet

Phosphotyrosine monoclonal antibody (orb1089153)

Description	Phosphotyrosine monoclonal antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Other, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated Phosphotyrosine
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.