

Product Datasheet

Strep-Tag II monoclonal antibody (orb1089166)

Description	Strep-Tag II monoclonal antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide Strep-Tag
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
Dilution Range	WB=1:1000-5000
Expiration Date	12 months from date of receipt.