

Product Datasheet

BRD1 Polyclonal Antibody (orb1089222)

Description	BRD1 Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Bovine, Canine, Equine, Gallus, Human, Mouse, Rabbit, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human BRD1:701-800/1058
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	140 kD
Uniprot ID	O95696
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.