

Product Datasheet

ERR β polyclonal antibody [KU] (orb1089281)

Description	ERR β polyclonal antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Synthetic peptide, corresponding to Human ERR β .
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 48 kDa
Uniprot ID	O95718
Dilution Range	WB: 1:1000~1:2000 IF: 1:50~1:200
Expiration Date	12 months from date of receipt.