


Product Datasheet

3-Nitrotyrosine antibody (orb108935)

Description	100 µg monoclonal antibody in 20 mM sodium phosphate, 150 mM sodium chloride, 50% glycerol, 3mM sod...
Species/Host	Mouse
Reactivity	All
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	3-Nitrotyrosine-KLH
Concentration	See vial for concentration
Storage	Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles
Note	For research use only
Application notes	Antibody can be used for Western blotting (Suggested dilution - 1:1000) and immunohistochemistry. Optimal concentration should be evaluated by serial dilutions.
Isotype	IgG1
Clonality	Monoclonal
Source	Mouse monoclonal antibody 2A12, obtained after immunization of mice with 3-Nitrotyrosine-KLH
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse brain lys...