

Product Datasheet

IL13 antibody (orb10895)

Description

IL13 antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IHC-P

Applications
Immunogen

KLH conjugated synthetic peptide derived from human IL13. Please contact us for the exact immunogen sequence. The peptide is available as orb374729.

Target

IL13

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

12 kDa

Uniprot ID

[P20109](#), [P35225](#), [P42203](#)

NCBI


[002188](#), [21](#), [002179](#)

Entrez
[3596](#)
Dilution Range


IHC-P: 1:50-300 (based on 0.5 mg/ml)

Expiration Date


12 months from date of receipt.



IHC-P image of rat epididymis tissue using...



Immunohistochemical staining of paraffin...



Immunohistochemical staining of paraffin...