

Product Datasheet

FOXP1 antibody (orb1089920)

Description	Human polyclonal antibody to FOXP1
Species/Host	Human
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IP
Immunogen	This antibody was obtained by recombinant antibody (rAb) phage display recognizing FOXP1 protein under non-denaturing conditions.
Target	FOXP1
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Application notes	This antibody has been used and recommended for the following methods: ELISA; spiked-IP; IF; FC.
Isotype	Fab fragment
Clonality	Monoclonal
Purity	Purified
Uniprot ID	Q9H334
Expiration Date	12 months from date of receipt.