

Product Datasheet

SETD7 antibody (orb1089925)

Description	Human polyclonal antibody to SETD7
Species/Host	Human
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ChIP, ELISA, FC, IF, IP
Immunogen	This antibody was obtained by recombinant antibody (rAb) phage display recognizing SETD7 protein under non-denaturing conditions.
Target	SETD7
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Application notes	This antibody has been used and recommended for the following methods: ELISA; IF; spiked IP; IP-MS; ChIP; FC.
Isotype	Fab fragment
Clonality	Monoclonal
Purity	Purified
Uniprot ID	Q8WTS6
Expiration Date	12 months from date of receipt.