

Product Datasheet

dsRNA antibody (orb1089973)

Description	Mouse monoclonal antibody to dsRNA
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	RDV-RNA-methylated bovine serum albumin complex.
Target	dsRNA
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Application notes	The antibody binds specifically to ds-RNA, which is part of the genetic material of some viruses, such as rotaviruses and the bluetongue virus. The antibody can be used as a diagnostic tool, to determine the nature of an unknown pathogen, and could potentially be used in gene therapy.
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.