


Product Datasheet

UEV1A antibody [Out of stock] (orb109002)

Description	Rabbit polyclonal to UEV1A
Species/Host	Rabbit
Tested Applications	WB
Immunogen	Rabbit polyclonal UEV1A antibody was raised against a peptide corresponding to 15 amino acids near the C-terminus of human UEV1A.
Preservatives	Provided as solution in phosphate buffered saline with 0.02% sodium azide
Note	For research use only
Isotype	IgG
Uniprot ID	Q13404
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat cell lys...