

Product Datasheet

Orexin Receptor antibody (orb109015)

Description	Rabbit polyclonal antibody to Orexin Receptor
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Synthetic peptide made to recombinant Human Orexin Receptor Type 2
Preservatives	Affinity purified polyclonal antibody provided as solution in phosphate buffered saline with 0.08% sodium azide
Form/Appearance	Affinity purified polyclonal antibody provided as solution in phosphate buffered saline with 0.08% sodium azide
Concentration	See vial for concentration
Storage	Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles
Note	For research use only
Application notes	Optimal concentration should be evaluated by serial dilutions.
Isotype	IgG
Uniprot ID	O43614
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	12 months from date of receipt.