


Product Datasheet

Recombinant PD-L1 antibody (orb1090509)


Description	Rabbit monoclonal antibody to PD-L1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Target	CD274
Form/Appearance	Purified recombination monoclonal antibody supplied in PBS with 0.05% (W/V) sodium azide, and 0.05% rAlbumin. This antibody is purified through a protein A column.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	33275
Uniprot ID	Q9NZQ7
Dilution Range	WB: 1:1000, IHC-P: 1:100-1:400
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...



Immunohistochemical analysis of paraffin...



Immunohistochemical analysis of paraffin...