



Product Datasheet

IL8 antibody (orb10915)


Description	IL8 antibody
Species/Host	Rabbit
Reactivity	Human, Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human IL-8 (31-99/99 aa)
Target	IL-8/CXCL8
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	8/9 kDa
Uniprot ID	P10145
Expiration Date	12 months from date of receipt.



Flow cytometric analysis of human hepatocellular carcinoma cells.



FACS of human hepatocellular carcinoma cells.



Immunohistochemical analysis of formalin-fixed paraffin-embedded tissue.