


Product Datasheet

Cytochrome P450 2A6 antibody (orb109150)

Description

Provided as solution in phosphate buffered saline



Western blot analysis of recombinant cyt...

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested
Applications

ELISA, IHC-P, WB

Concentration

See vial for concentration

Storage


Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles

Note

For research use only

Application notes

Antibody can be used for Western blotting (1-2 µg/ml), immunofluorescence and immunohistochemistry (formalin-fixed, paraffin-embedded tissues). ELISA titer to be determined by end-user. Optimal concentration should be evaluated by serial dilutions.



Immunohistochemical staining of human li...

Isotype

IgG1

Clonality

Monoclonal

Source

Hybridoma produced by the fusion of splenocytes from BALB/c mice immunized with Synthetic peptide conjugated to ovalbumin derived from the C-terminus of CYP2A6 and mouse myeloma Ag8563 cells. Antigen identical in CYP2A6, 2A7 and 2A13.

Uniprot ID

P05177

Hazard
Information

This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.

Expiration Date

12 months from date of receipt.