



Product Datasheet

HDAC2 antibody (orb10916)

Description	Rabbit polyclonal antibody to HDAC2
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human HDAC2(NP_001518.3).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	54kDa
Uniprot ID	Q92769
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200, IF/ICC: 1:50-1:100, IP: 1:20-1:50, FACS: 1:20-1:50
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human lung tissue.



Western blot analysis of extracts of H460 cells.