

Product Datasheet

Serine Palmitoyltransferase 1 antibody (orb109166)

Description

Rabbit polyclonal to Serine


Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC, WB

Applications

Immunohistochemical staining of human lu...

Immunogen


Synthetic peptide derived from the human serine palmitoyltransferase enzyme.

Preservatives

Provided as ligand affinity purified antibody in phosphate buffered saline with 0.08% sodium azide

Form/Appearance

Provided as ligand affinity purified antibody in phosphate buffered saline with 0.08% sodium azide



Western blot analysis of human kidney ly...

Concentration

See vial for concentration

Storage

Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles

Note

For research use only

Application notes

Antibody can be used for Western blotting (1-5 µg/ml) and immunohistochemistry (5-10 µg/ml). Optimal concentration should be evaluated by serial dilutions.

Isotype

IgG

Uniprot ID
O15269
Hazard
Information

This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.

Expiration Date

12 months from date of receipt.