



Product Datasheet

Insulin antibody (orb10922)


Description	Insulin antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human Insulin (46-59/110 aa)
Target	Proinsulin
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	5.8/12 kDa
CAS Number	11061-68-0
Uniprot ID	P01308
Dilution Range	IHC-P=1:100-500, IF=1:100-500, IHC-F=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



IHC-p on mouse pancreatic glands with or...



Insulin expression in mouse pancreatic g...



Immunohistochemical analysis of formali...