



Product Datasheet

RPA1 antibody (orb10926)

Description	Rabbit polyclonal antibody to RPA1
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Recombinant fusion protein of human RPA70/RPA70/RPA1(NP_002936.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	68kDa
Uniprot ID	P27694
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200, IF/ICC: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of MCF-7 cells...



Western blot analysis of extracts of various tissues...