


Product Datasheet

IRS4 antibody (orb10939)

Description	IRS4 antibody
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	IHC-Fr
Immunogen	KLH conjugated synthetic peptide derived from mouse IRS-4 (401-500/1216 aa)
Target	IRS-4
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	134 kDa
Uniprot ID	Q9Z0Y7
Dilution Range	IHC-F=1:100-500, IHC-P=1:100-500, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



IHC-P of rat pancreatic tissue using IRS...