

## Product Datasheet

### Phospho-GSK3 beta (Ser25) Antibody (orb1093956)

|                                |                                       |
|--------------------------------|---------------------------------------|
| <b>Description</b>             | Phospho-GSK3 beta<br>(Ser25) Antibody |
| <b>Species/Host</b>            | Rabbit                                |
| <b>Reactivity</b>              | Human, Mouse, Rat                     |
| <b>Conjugation</b>             | Unconjugated                          |
| <b>Tested<br/>Applications</b> | ELISA, IHC-P, WB                      |
| <b>Note</b>                    | For research use only                 |
| <b>Clonality</b>               | Polyclonal                            |
| <b>Expiration Date</b>         | 12 months from date of<br>receipt.    |