

Product Datasheet

Sphingomyelin Phosphodiesterase 3 antibody (orb109408)

Description

Rabbit polyclonal to Sphingomyelin

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

WB

Applications

Immunogen

Synthetic peptide derived from human sphingomyelin phosphodiesterase 3

Preservatives

Provided as solution in phosphate buffered saline with 0.08% sodium azide

Form/Appearance

Provided as solution in phosphate buffered saline with 0.08% sodium azide

Concentration

Lot specific, see vial

Storage

Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles

Note

For research use only

Application notes

Antibody can be used for Western blotting (1:400 starting dilution). Optimal concentration should be evaluated by serial dilutions. Protein appears to run lower than expected protein molecular weight. However, antibody is specific for the protein and band blocks when blocking peptide used.

Uniprot ID

Q9NY59


Hazard

Information

This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.

Expiration Date

12 months from date of receipt.



Western blot analysis of human brain lys...