


Product Datasheet

PCSK9 antibody (orb109430)

Description	Rabbit polyclonal to PCSK9
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Synthetic peptide derived from the human PCSK9 protein reported to block PCSK9-LDLR interaction.
Preservatives	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Form/Appearance	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Concentration	See vial for concentration
Storage	Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles
Note	For research use only
Application notes	Optimal concentration should be evaluated by serial dilutions.
Isotype	IgG
Uniprot ID	Q8NBP7
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Hep G2 C...