



Product Datasheet

Phospho-TFEB (Ser211) Antibody (orb1095423)


Description	Phospho-TFEB (Ser211) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Human kidney tissue



Immunofluorescence image of A549 cells



Immunohistochemical staining of Human kidney tissue