

## Product Datasheet

### Phospho-AQP2 (Ser256) Antibody (orb1095458)

|                            |                                   |
|----------------------------|-----------------------------------|
| <b>Description</b>         | Phospho-AQP2 (Ser256)<br>Antibody |
| <b>Species/Host</b>        | Rabbit                            |
| <b>Reactivity</b>          | Human, Mouse, Rat                 |
| <b>Conjugation</b>         | Unconjugated                      |
| <b>Tested Applications</b> | ELISA                             |
| <b>Note</b>                | For research use only             |
| <b>Clonality</b>           | Polyclonal                        |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |