

Product Datasheet

Small ubiquitin-related modifier 1 antibody (orb109551)

Description

Rabbit polyclonal to Small ubiquitin-related


Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunohistochemical staining of normal h...

Immunogen

Synthetic peptide derived from the human SUMO-1 protein.

Preservatives

Provided as solution in phosphate buffered saline with 0.08% sodium azide

Form/Appearance

Provided as solution in phosphate buffered saline with 0.08% sodium azide

Concentration

See vial for concentration

Storage

Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles

Note

For research use only

Application notes

Antibody can be used for Western blotting (1:400 dilution) and immunohistochemistry (10-15µg/ml). Optimal concentration should be evaluated by serial dilutions.

Isotype

IgG

Uniprot ID
P63165
Hazard
Information

This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.

Expiration Date

12 months from date of receipt.