

## Product Datasheet

### Goat Anti-Human IgG, Fcγ fragment specific Antibody (orb1095814)

|                            |                                                               |
|----------------------------|---------------------------------------------------------------|
| <b>Description</b>         | Goat Anti-Human IgG, Fcγ fragment specific Antibody           |
| <b>Species/Host</b>        | Goat                                                          |
| <b>Reactivity</b>          | Human                                                         |
| <b>Conjugation</b>         | Unconjugated                                                  |
| <b>Tested Applications</b> | ELISA                                                         |
| <b>Immunogen</b>           | IgG, Fcγ fragment specific                                    |
| <b>Preservatives</b>       | 0.01M Sodium Phosphate,<br>0.25M NaCl, pH 7.6                 |
| <b>Form/Appearance</b>     | Liquid                                                        |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| <b>Note</b>                | For research use only                                         |
| <b>Isotype</b>             | IgG                                                           |
| <b>Clonality</b>           | Polyclonal                                                    |
| <b>Expiration Date</b>     | 12 months from date of receipt.                               |