

Product Datasheet

N Antibody (orb1095831)

Description	N Antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human Novel Coronavirus Nucleoprotein (N) (1-419aa)
Target	N
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P0DTC9
Dilution Range	WB:1000-5000, ELISA:1:1000-1:5000
Expiration Date	12 months from date of receipt.