

Product Datasheet

TNFRSF17 Antibody (orb1095850)

Description	TNFRSF17 Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human TNFRSF17 protein (1-54AA)
Target	TNFRSF17
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Uniprot ID	Q02223
Dilution Range	WB: 1:1000-1:5000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.