

Product Datasheet

TSG101 Antibody (orb1095854)

Description	TSG101 Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, WB
Immunogen	Recombinant Human TSG101 protein (1-145AA)
Target	TSG101
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	Q99816
Dilution Range	WB: 1:500-1:2000, IF: 1:50-1:200, FC: 1:50-1:200
Expiration Date	12 months from date of receipt.