

Product Datasheet

GST antibody (orb109587)

Description	Mouse monoclonal to GST
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	A fusion protein to the internal sequence GST
Preservatives	Provided as 0.2µm sterile filtered solution in phosphate buffered saline with 0.08% sodium azide.
Form/Appearance	Provided as 0.2µm sterile filtered solution in phosphate buffered saline with 0.08% sodium azide.
Concentration	See vial for concentration
Storage	Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles
Note	For research use only
Application notes	Application: Western Blots at 1:1000- 1:2,000 Optimal concentration should be evaluated by serial dilutions.
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P10649
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	12 months from date of receipt.