

Product Datasheet

Recombinant Human Osteopontin(SPP1),partial (orb1096381)

Description	Recombinant Human Osteopontin(SPP1),partial
Target	SPP1
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-tagged
Note	For research use only
Protein Sequence	IPVKQADSGSSEEKQLYNKYPDAVATWLNPDPSQKQNLAPQNAVSSSEETNDFK QETLPSKSNESHDMDDMDEDDDDHVDSQDSIDSNDSDDVDDTDDSHQSD ESHHSDESDELVTDFPTDLPATEVFTPVPTVDYDGRGDSVVYGLRSKSKKFR RPDIQYPDATDEDITSHMESEELNGAYKAIPVAQDLNAPSDWDSRGKDSYETSQ LDDQSAETHSHKQSRLYKRKANDESNEHSDVIDSQELSKVSREFHSHEFHSHED MLVVDPKSKEEDKHLKF
Purity	Greater than 90% as determined by SDS-PAGE.
Source	Yeast
MW	34.2 kDa
Uniprot ID	P10451
Expiration Date	6 months from date of receipt.