


Product Datasheet

LDL Receptor antibody (orb10972)


Description	LDL Receptor antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, ICC, WB
Immunogen	KLH conjugated synthetic peptide derived from human LDL-R (781-860/860 aa)
Target	LDL receptor
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	92 kDa
Uniprot ID	P01130
Dilution Range	WB=1:500-2000, ICC=1:100, Flow-Cyt=1µg/Test
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...



Western blot of rat brain tissue lysate ...



IHC-P of rat brain tissue (orb10972 at 1...