

Product Datasheet

p90RSK (phospho-Thr359/Ser363) rabbit pAb Antibody (orb1097619)

Description	p90RSK (phospho-Thr359/Ser363) rabbit pAb Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q15418
Expiration Date	12 months from date of receipt.