

## Product Datasheet

# Ribosomal Protein S6 (phospho Ser240) Polyclonal Antibody (orb1097697)

|                            |                                                                                              |
|----------------------------|----------------------------------------------------------------------------------------------|
| <b>Description</b>         | Ribosomal Protein S6 (phospho Ser240) Polyclonal Antibody                                    |
| <b>Species/Host</b>        | Rabbit                                                                                       |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                            |
| <b>Conjugation</b>         | Unconjugated                                                                                 |
| <b>Tested Applications</b> | ELISA, IF, IHC-P, WB                                                                         |
| <b>Storage</b>             | Store at -20°C for long term preservation (recommended).<br>Store at 4°C for short term use. |
| <b>Note</b>                | For research use only                                                                        |
| <b>Clonality</b>           | Recombinant                                                                                  |
| <b>Uniprot ID</b>          | <b>P62753</b>                                                                                |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                              |