

## Product Datasheet

### Aquaporin 4 Monoclonal Antibody [Out of stock] (orb1097862)

|                            |                                 |
|----------------------------|---------------------------------|
| <b>Description</b>         | Aquaporin 4 Monoclonal Antibody |
| <b>Species/Host</b>        | Mouse                           |
| <b>Tested Applications</b> | IF, WB                          |
| <b>Storage</b>             | Store at -20°C/1 year           |
| <b>Note</b>                | For research use only           |
| <b>Clonality</b>           | Monoclonal                      |
| <b>Uniprot ID</b>          | <b>P55087</b>                   |
| <b>Expiration Date</b>     | 12 months from date of receipt. |