

Product Datasheet

Leptin Receptor antibody (orb10980)

Description

Leptin Receptor antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Porcine, Rat

Conjugation

Unconjugated

Tested

ICC, IF, IHC-P, WB

Applications

Immunogen

KLH conjugated synthetic peptide derived from human Leptin Receptor. Please contact us for the exact immunogen sequence. The peptide is available as orb374731.

Target

Leptin Receptor

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

132 kDa

Uniprot ID

[P48356](#), [Q62959](#), [P48357](#)

NCBI

[31](#), [001003679](#), [11](#)

Entrez


[3953](#)

Dilution Range


WB: 1:100-1000, IHC-P: 1:100-500, IF/ICC: 1:100-800 (based on 0.5 mg/ml)

Expiration Date


12 months from date of receipt.



Immunofluorescence analysis of pig small...



IHC-P staining of pig small intestines t...



IF image of pig small intestines tissue ...