

## Product Datasheet

### **delta 1 Catenin/CAS/CTNND1 Antibody (orb1098013)**

## Description

delta 1 Catenin/CAS/CTNND1 Antibody

### Species/Host

Rabbit

### Reactivity

Human, Mouse, Rat

### Conjugation

Unconjugated

### Tested

ELISA, FC, WB

### Applications

### Immunogen

E.coli-derived human delta 1 Catenin/CAS/CTNND1 recombinant protein (Position: E31-R875).

### Form/Appearance

Lyophilized

### Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

### Storage

At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.

### Note

For research use only

### Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

### Isotype

Rabbit IgG

### Clonality

Polyclonal

### MW

100-110 kDa

### Uniprot ID


[O60716](#)

### Dilution Range


Western blot, 0.25-0.5 µg/ml, Human, Mouse, Rat  
Flow Cytometry, 1-3 µg/1x10<sup>6</sup> cells, Human Direct  
ELISA, 0.1-0.5 µg/ml, Human

### Expiration Date

12 months from date of receipt.



Western blot analysis of Delta 1 Catenin...



Flow Cytometry analysis of Caco-2 cells ...