


Product Datasheet

Human Leptin ELISA Kit (DIY Antibody Pairs) (orb1098131)

Description	Human Leptin ELISA Kit (DIY Antibody Pairs)
Reactivity	Human
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: V22-C167
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	18641 MW
Uniprot ID	A4D0Y8
Expiration Date	Please enquire.



Human
Leptin EZ
Set ELISA
Kit
standard
C...