


Product Datasheet

Human IL-5 ELISA Kit (DIY Antibody Pairs) (orb1098137)

Description	Human IL-5 ELISA Kit (DIY Antibody Pairs)
Reactivity	Human
Range	7.8 pg/ml - 500 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: sf21; Immunogen sequence: I20-S134
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) and urine.
MW	15238 MW
Uniprot ID	P05113
Expiration Date	Please enquire.



Human IL-5
EZ Set
ELISA Kit
standard
cur...