


Product Datasheet

Human IL-32 ELISA Kit (DIY Antibody Pairs) (orb1098148)

Description	Human IL-32 ELISA Kit (DIY Antibody Pairs)
Reactivity	Human
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: A31-K234
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA).
Uniprot ID	P24001
Expiration Date	Please enquire.



Human IL-32 EZ-Set ELISA Kit standard cu...