

Product Datasheet

Human Cystatin C ELISA Kit (DIY Antibody Pairs) (orb1098151)

Description Human Cystatin C ELISA Kit (DIY Antibody Pairs)

Reactivity Human

Range 312 pg/ml - 20,000 pg/ml

Tested Applications ELISA

Immunogen Expression system for standard: NS0;
Immunogen sequence: M1-A146


Storage Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)

Note For research use only

Sample Types cell culture supernatants, serum, plasma (heparin, EDTA) , saliva, urine and human milk.

Uniprot ID **P01034**

Expiration Date Please enquire.



Human
Cystatin C
EZ-Set
ELISA Kit
standa...