


Product Datasheet

Human Trail ELISA Kit (DIY Antibody Pairs) (orb1098163)

Description	Human Trail ELISA Kit (DIY Antibody Pairs)
Reactivity	Human
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: T95-G281
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum, plasma (heparin, EDTA) and saliva.
Uniprot ID	P50591
Expiration Date	Please enquire.



Human
Trail EZ-Set
ELISA Kit
standard
cu...