


Product Datasheet

Mouse TFF1 ELISA Kit (DIY Antibody Pairs) (orb1098165)

Description	Mouse TFF1 ELISA Kit (DIY Antibody Pairs)
Reactivity	Mouse
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: Q22-F87
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
Uniprot ID	Q08423
Expiration Date	Please enquire.



Mouse
TFF1 EZ-
Set ELISA
Kit
standard
cur...