


Product Datasheet

Human IGFBP3 ELISA Kit (DIY Antibody Pairs) (orb1098166)

Description	Human IGFBP3 ELISA Kit (DIY Antibody Pairs)
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: A28Q-291
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
Uniprot ID	P17936
Expiration Date	Please enquire.



Human
IGFBP3 EZ-
Set ELISA
Kit
standard
C...