


Product Datasheet

Human FABP4 ELISA Kit (DIY Antibody Pairs) (orb1098170)

Description	Human FABP4 ELISA Kit (DIY Antibody Pairs)
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: M1-A132
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) and cell lysates.
Uniprot ID	P15090
Expiration Date	Please enquire.



Human
FABP4 EZ-
Set ELISA
Kit
standard
CU...