


Product Datasheet

Human GDF-15 ELISA Kit (DIY Antibody Pairs) (orb1098177)

Description	Human GDF-15 ELISA Kit (DIY Antibody Pairs)
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: CHO; Immunogen sequence: A197-I308
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
Uniprot ID	Q99988
Expiration Date	Please enquire.



Human
GDF-15 EZ-
Set ELISA
Kit
standard
C...