


Product Datasheet

Human TREML1 ELISA Kit (DIY Antibody Pairs) (orb1098178)

Description	Human TREML1 ELISA Kit (DIY Antibody Pairs)
Reactivity	Human
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Q16-P16
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
Uniprot ID	Q86YW5
Expiration Date	Please enquire.



Human
TREML1
EZ-Set
ELISA Kit
standard
C...